

TECHNICAL DATA SUMMARY (TDS)

VACLOK SYRINGES

1. DEVICE DESCRIPTION

The VacLok® syringe is a modification to Merit's Medallion syringe family. The VacLok® syringe creates a negative pressure or vacuum for applications that are typically performed in the cardiac catheterization and radiology labs. The 10, 20, 30 and 60 ml VacLok® syringes are designed to stop in multiple positions when pulling a vacuum and have the capability of holding a negative pressure by engaging a stop pin that fits into the syringe barrel.

The VacLok® syringe barrels are manufactured of clear polycarbonate and are available with a fixed male Luer lock connector. VacLok barrels are created to house a stop pin that can lock the plunger at set intervals. The plunger and stop pin are molded from Acrylonitrile-Butadiene-Styrene (ABS) and are fitted with a silicone rubber piston seal. The 60 mL VacLok® syringe utilizes an EPDM O-Ring seal instead of the silicone piston seal used in other syringes.



2. INTENDED USE / INDICATION FOR USE

The Merit Syringe is intended to be used by a cardiologist or radiologist during angiographic or radiologic procedures.

3. Contraindications

There are no contraindications or warnings for this product.

4. INSTRUCTIONS FOR USE

The Instructions for Use and product labelling for the product have been formatted in multiple languages.

5. Classification:

a. European Community:

The classifications take into consideration of all rules based on the component with the highest classification. For devices claiming conformity to the Essential Requirements of the MDR, the following classification rules are used:

The product is not an active device, and none of the special rules apply.

| Subjects | Applied Rules | Justification |
|----------------------|--|----------------------------------|
| Non-invasive devices | <input type="checkbox"/> 1, <input checked="" type="checkbox"/> 2, <input type="checkbox"/> 3, <input type="checkbox"/> 4 | See table below. |
| Invasive devices | <input type="checkbox"/> 5, <input type="checkbox"/> 6, <input type="checkbox"/> 7, <input type="checkbox"/> 8 | Not an invasive device. |
| Active devices | <input type="checkbox"/> 9, <input type="checkbox"/> 10, <input type="checkbox"/> 11, <input type="checkbox"/> 12 | Not an active device. |
| Special rules | <input type="checkbox"/> 13, <input type="checkbox"/> 14, <input type="checkbox"/> 15, <input type="checkbox"/> 16, <input type="checkbox"/> 17, <input type="checkbox"/> 18 | None of the special rules apply. |

Table 1: MDR Classification Rationale of VacLok

| Component | Rationale | Rule | Class |
|-----------|--|------|-----------------------|
| Syringe | Non-invasive devices intended for channeling liquids... for the purpose of eventual infusion, administration or introduction into the body | 2 | I – Sterile Measuring |

b. United States Food and Drug Administration:

Class: II
 Product Code/CFR: FMF / 21 CFR 880.5860
 510(k) Reference: K994253
 Device Name: Piston Syringe

c. Canada:

Class: II (Rule 5)
 License: 6070

d. GMDN and CND

GMDN Code: 47017 - General-purpose syringe, Single Use
 CND Code: A0201 - Single Use Syringes

6. DEVICE SPECIFICATIONS

| SPECIFICATION NAME | DESCRIPTION |
|--------------------|--|
| Capacity | 10mL, 20mL, 30mL and 60mL |
| Luer Connections | Fixed Male Luer & Male Slip Luer Per ISO 594-1/2 (or ISO 80369-7 testing currently underway) |
| Graduation Scale | Per ISO 7886-1 |
| Locking mechanism | Four position: 10ml and 20mL Six position: 30mL and 60mL |

7. MATERIALS OF CONSTRUCTION

| PRODUCT COMPONENT | MATERIAL |
|-------------------|---|
| Barrel | Polycarbonate |
| Stop Pin | Acrylonitrile-Butadiene-Styrene (ABS) |
| Plunger | Acrylonitrile-Butadiene-Styrene (ABS) |
| 10-30 mL Tip | Silicone rubber |
| O-ring (60mL) | Ethylene Propylene Diene Methylene (EPDM) |
| Ink | Black or White U.V. Curable |
| Lubricant | Medical grade silicone fluid |

| PACKAGING COMPONENT | MATERIAL |
|---------------------|--|
| Blister | Formed blister of 1073 Tyvek and EVA Trilayer film |
| Carton | Corrugate single wall |

8. Packaging

Merit's syringes are packaged as sterile, single use devices. The barrier properties of the packaging assure that sterility is maintained. All units are packaged in a pre-formed blister tray and sealed with medical grade Tyvek™, or a Nylon pouch configuration. VacLok syringes are packaged 25 units per carton for 10 mL & 60 mL and 20 units per carton for 20 mL & 30 mL, with 4 cartons per shipper.

Non-sterile syringe assemblies are packaged to facilitate incorporation into customized procedure packs which are provided to clinicians as sterile, single use devices. Product labeling of non-sterile devices clearly specifies the requirement for further processing (e.g. sterilization).


Certain bulk non-sterile syringe assemblies include labeling that specifies the validated sterilization parameters that repackers can use to leverage Merit's CE mark. These are triple poly-bagged and placed into cartons to protect the product and to reduce particulate contamination

Figure 1: Syringe Line Multivac Packaging



9. Labeling

Figure 2: Unit Label - VacLok



VacLok® Syringe

REF VAC120
VERSION A

LOT VAC120/A

Manufacture Date:
2022-05-18

Use By:
2022-05-18

STERILE EO



Single Use



Do not resterilize

Single sterile barrier system

Rx ONLY



Caution

Non-pyrogenic



Do not use if package is damaged and consult instructions for use

MD Medical Device

Made in U.S.A. **1 Unit**

For Patent Coverage, See www.merit.com/patents



Merit Medical Systems, Inc., 1600 W. Merit Parkway, South Jordan, Utah 84095
U.S.A. Customer Service 1-800-356-3748

EC REP Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
European Customer Service +31 43 358 82 22

i Consult Instructions for Use. For electronic copy scan QR code,
or go to www.merit.com/ifu and enter IFU ID:1545494
For printed copy, call U.S.A. or EU Customer Service



=CL(0:0:5:0)<HH=MS(9:0)<MO/DD/YYYY>CN(0:0:5:+1:3)0000FSYM0016 REV 001 02/21



Figure 3: Carton Label - VacLok

REF VAC120
VERSION A

LOT VAC120/A

Manufacture Date:
2022-05-18

Use By:
2022-05-18

Single Use

Do not use if package is damaged and consult instructions for use


STERILE EO Non-pyrogenic For Patent Coverage, See www.merit.com/patents

Caution Do not resterilize **MD** Medical Device **Rx ONLY** **20 Units**

Merit Medical Systems, Inc., 1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
1-801-253-1600 U.S.A. Customer Service 1-800-356-3748


EC REP Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
European Customer Service +31 43 358 82 22

Consult Instructions for Use. For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID:1545494
For printed copy, call U.S.A. or EU Customer Service




20ml ± 4%

VacLok® Syringe



000000000001

UDI



(01)10884450610654

(01)10884450610654(17)220518(10)VAC120/A

ID 080114 403135001/A

MERITMEDICAL®

Made in U.S.A.

CE 2797

001 05/18/22 HH:mm:ss LFSYZ0013 REV 001 02/21

Figure 4: Shipper Label - VacLok

REF VAC120
VERSION A

LOT VAC120/A

Manufacture Date:
2022-05-18

Use By:
2022-05-18

20ml ± 4%

VacLok® Syringe

Single Use

Do not use if package is damaged and consult instructions for use

STERILE EO Non-pyrogenic For Patent Coverage, See www.merit.com/patents

Caution Do not sterilize **MD** Medical Device **Rx ONLY** 80 Units

Merit Medical Systems, Inc., 1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
European Customer Service +31 43 358 82 22

Consult Instructions for Use. For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID:1545494
For printed copy, call U.S.A. or EU Customer Service

Made in U.S.A.
CE 2797

UDI
0000000000001
0000000000001
(01)20884450610651
(01)20884450610651(17)220518(10)VAC120/A

MERITMEDICAL®

001 05/18/22 HH:mm:ss LFSYZ0013 REV 001 02/21

10. Compliance with Relevant Regulations and Technical Standards

EU:

Medical Device Regulation (MDR) of the European Union (Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices). The applicable Certificate is MDR 716389.

US:

The product is legally marketed according to the requirements of the United States Food and Drug Administration in compliance with the premarket notification 510(k) process. The 510(k) clearance letter is K994253.

CAN:

Class: II (Rule 5)
License: 6070

The qualification activities that were performed for the product were based on the below applicable standards:

| DOCUMENT | TITLE |
|---------------------------------------|---|
| Horizontal (General) Standards | |
| Council Directive 93/42/EEC | Medical Device Directive of the European Union |
| Regulation (EU) 2017/745 | Medical Device Regulation (MDR) of the European Union (Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices) |
| ISO 13485 | Quality Systems – Medical Devices – Quality Management Systems. Requirements for Regulatory Purposes |
| MEDDEV 2.1/5 (June 1998) | Medical Devices: Guidance Document: Medical Devices with a Measurement Function |
| ISO 14971 | Medical devices — Application of risk management to medical devices |
| ISO 10993-1 | Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process |
| ISO 11607-1 | Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems, and packaging systems. |
| ISO 11607-2 | Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes |
| ASTM D4169 | Standard practice for performance testing of shipping containers and systems |
| ASTM F1980 | Standard guide for accelerated aging of sterile barrier systems for medical devices |
| ASTM F2475 | Standard Guide for Biocompatibility of Medical Device Packaging Materials |
| ASTM F 2096 | Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurization (Bubble Test) |
| ASTM F 1929 | Standard Test Method for Detecting Seal Leaks in porous Medical Packaging by Dye Penetration |
| ASTM F88/F88M | Standard Test Method for Seal Strength of Flexible Barrier Materials. |
| BS EN 980 | Symbols for use in the labeling of medical devices |
| BS EN 1041 | Information supplied by the manufacturer of medical devices |
| ISO 15223-1 | Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements |

| DOCUMENT | TITLE |
|-----------------------------------|--|
| ISO 11135-1 | Sterilization of health care products -- Ethylene oxide -- Requirements for development, validation and routine control of a sterilization process for medical devices |
| EN556-1 | Sterilization of medical devices – Requirements for terminally-sterilized medical devices to be labeled “Sterile” |
| USP | United States Pharmacopeias |
| AAMI TIR 28 | Product Adoption and process equivalency for ethylene oxide sterilization |
| ANSI/AAMI ST72:2011 | Bacterial Endotoxins – Test methods, routine monitoring, and alternatives to batch testing |
| Product Specific Standards | |
| ISO 7886-1 | Sterile hypodermic syringes for single use – Part 1: Syringes for manual use |
| ISO 594-1 | Conical Fittings with 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment – Part 1: General Requirements |
| ISO 594-2 | Conical Fittings with 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment – Part 2: Lock Fittings |
| ISO 80369-7 | Small bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications |

11. Storage and handling

Store under general warehouse conditions.

12. Biocompatibility

The Product meets ISO 10993-1 and the FDA guidance document: Use of International Standards ISO-10993-1, “Biological Evaluation of Medical Devices Part 1: Evaluation and Testing Within a Risk Management Process.”

13. SHELF LIFE

The Products are labelled for a 5-year shelf life.
Non-sterile units are labeled for 5-year shelf life.

14. STERILIZATION DETAILS

Merit Medical Systems, Inc. utilizes Ethylene Oxide (EO) sterilization for the products described in this Technical Data Summary. The sterilization cycle is validated according to European International Standard ISO 11135-1. These products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals.

15. MATERIAL COMPLIANCE

EU REACH: To the best of our knowledge based on information currently available from our raw material suppliers as well as an assessment of the general material properties, the product does not contain substances listed on the ECHA Candidate List of Substances of Very High Concern (SVHC).

CMR and Endocrine-Disrupting Substances: An assessment of the raw materials used to manufacture the product in this summary with respect to carcinogenic, mutagenic or toxic to reproductive ('CMR') compounds, of category 1A or 1B, in accordance with Part 3 of Annex VI to Regulation (EC) No 1272 | 2008 of the European Parliament and of the Council, or substances having endocrine-disrupting properties was completed. No CMR compounds or substances having endocrine-disrupting properties are contained in the products that are the subject of this Technical Data Summary.

16. MANUFACTURING CONDITIONS / QUALITY CONTROL

Merit Medical Systems, Inc. manufactures the products in this summary in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820 *Quality Systems Regulations*.

Merit Medical Systems facilities used to manufacture the product in this summary are also certified to ISO 13485 Quality Management System.