

Guide Wire Accessories



BowTie™ Guide Wire
Insertion Tool Fits Guide
Wires ≤ 0.018" (0.46mm)

502423001

Metal Guide Wire Insertion Tool
Hub Material: Acrylic
Inner Lumen: 0.022"
Shaft Material: Stainless Steel

501994001

Plastic Guide Wire Insertion Tool
Hub Material: High Density Polyethylene
Inner Lumen: 0.028"
Shaft Material: High Density Polyethylene

MAP600B

Peel-Away Guide
Wire Insertion Tool,
Tapered Tip,
Inner Lumen 0.043"

502199002

Peel-Away Guide
Wire Insertion Tool,
Non-tapered Tip,
Inner Lumen 0.055"

502199003



Natural Color Retention Clip

100730002

J-Straightener
for .035" wires

101415001

J-Straightener
for .038" wires

101415002

Ringmaster™ Guide Wire Basin

RM01

2,500cc Guide Wire Bowl

M-P2000