

Adapters & Rotaters

90 Degree Right Angle

Material: Polycarbonate



Clear. Pocket ID .231"-.202"
Taper

100723001



Clear. Pocket ID .172"-.152"
Taper

100723002

Male Rotating

Body Material: Polycarbonate
O-Ring Material: EPDM



Clear. Tubing Pocket ID
.180"-.205"

510104002



Gamma. Tubing Pocket ID
.136"-.139"

501178001

Male/Female and Male/Male Rotating

Body Material: Polycarbonate
O-Ring Material: EPDM



Clear. Male/Female Rotating

562013001



Clear. Male/Male Rotating

MMA100P