

# CT Transfer Sets

Many tubing sets available with non-DEHP material.

Tubing Material: PVC

Inches: 20" O.D.: .160"



CT Contrast Transfer Set with  
Filtered Sosa Spike  
Swabbable Valve, ID: .120"

**TS120**



CT Large Bore Contrast Transfer Set  
Swabbable Valve, ID: .154"

**TS520L**



CT Saline Transfer Set with Blue Stripe Tubing and  
Longer Spike, Swabbable Valve, ID: .120"

**TS220**