

Contrast Transfer Sets

Many tubing sets available with non-DEHP material

Tubing Material: PVC
Inches: 20" O.D.: .160"



CT Contrast Transfer Set with
Filtered Sosa Spike
Swabbable Valve, ID: .120"

TS120



CT Large Bore Contrast Transfer Set
Swabbable Valve, ID: .154"

TS520L



CT Saline Transfer Set with Blue Stripe
Tubing and Longer Spike, Swabbable
Valve, ID: .120"

TS220



20 Small Bore, Female and Slip Luer

K08-01859P



Inline Contrast Saver, 54 Small Bore Tubing,
(2) 1-Way Check Valves, Std. Pressure
3-Way Stopcock, Male Luers

K08-00422D