

Coaxial Introducers



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK401



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK401-E



French Size: 4F
 Needle Length: 4cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK401-PED



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm
 Gauge: 21, Guide Wire Diameter/Length:
 0.018" (0.46mm) x 40cm

MAK401N



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm
 Gauge: 21, Guide Wire Diameter/Length:
 0.018" (0.46mm) x 40cm, Echo Enhanced

MAK401N-E



French Size: 4F
 Needle Length: 4cm
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm
 Gauge: 21, Guide Wire Diameter/Length:
 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK401N-PED

PERCUTANEOUS
ACCESS

Coaxial Introducers



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Platinum Tip
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm

MAK405



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Platinum Tip
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK405-E



French Size: 4F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK405-PED



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21,
 Guide Wire Diameter/Length: 0.018"
 (0.46mm) x 40cm

MAK501



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK501-E



French Size: 5F
 Needle Length: 4cm
 Core/Tip: Stainless Steel/Stainless Steel
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK501-PED

Coaxial Introducers



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Nitinol/Platinum Tip
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm

MAK501N



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Nitinol/Platinum Tip
 Length: 15cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm

MAK501N-15



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Nitinol/Platinum Tip
 Length: 10cm, Gauge: 21, Guide Wire
 Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK501N-E

PERCUTANEOUS ACCESS



French Size: 5F
 Needle Length: 4cm
 Core/Tip: Stainless Steel/Stainless Steel, Length: 10cm, Gauge: 21
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK501N-PED



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm

MAK505



French Size: 5F
 Needle Length: 7cm
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK505-E



French Size: 5F
 Needle Length: 4cm
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm
 Echo Enhanced

MAK505-PED