


# Drainage Catheters



## Why Use ReResolve® Locking Drainage Catheters?

- Highly flexible radiopaque marker band
- 25cm working length
- 0.038" (0.97mm) maximum guide wire
- Hydrophilic coating (20cm)
- Extra large drain holes

## Locking Drainage Catheters without Marker Band



6.5F, 7 Drainage Holes  
Green Hub

**RLC-6-038**



8.5F, 8 Drainage Holes  
Blue Hub

**RLC-8-038**



10F, 8 Drainage Holes  
Purple Hub

**RLC-10-038**



12F, 8 Drainage Holes  
Yellow Hub

**RLC-12-038**



14F, 9 Drainage Holes  
Orange Hub

**RLC-14-038**

## Locking Drainage Catheters with Marker Band



6.5F, 7 Drainage Holes  
Green Hub

**RLC-6-038MB**



8.5F, 8 Drainage Holes  
Blue Hub

**RLC-8-038MB**



10F, 8 Drainage Holes  
Purple Hub

**RLC-10-038MB**



12F, 8 Drainage Holes  
Yellow Hub

**RLC-12-038MB**



14F, 9 Drainage Holes  
Orange Hub

**RLC-14-038MB**

## Non-Locking Drainage Catheters with J Tip



J Tip, 6.5F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage  
Green Hub

**RNL-6-038J**



J Tip, 8.5F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage  
Blue Hub

**RNL-8-038J**



J Tip, 10F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage  
Purple Hub

**RNL-10-038J**



J Tip, 12F, 5 Drainage Holes  
Non-Locking (NL), Universal Drainage  
Yellow Hub

**RNL-12-038J**

# Non-Locking Drainage Catheters with Straight Tip



Straight Tip, 6.5F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage,  
Green Hub

**RNL-6-038S**



Straight Tip, 8.5F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage,  
Blue Hub

**RNL-8-038S**



Straight Tip, 10F, 4 Drainage Holes  
Non-Locking (NL), Universal Drainage,  
Purple Hub

**RNL-10-038S**



Straight Tip, 12F, 5 Drainage Holes  
Non-Locking (NL),  
Universal Drainage, Yellow Hub

**RNL-12-038S**