

Drip Chambers

Many tubing sets available with non-DEHP material.

Cap: Polyethylene Spike/Inline Tubing Adapter: ABS

Chamber: Polypropylene Vent: Polyethelyne



Macro Inline Non-Vented
700098



Drip Chamber, Vented with
15 Micron Filter
700026



Drip Chamber Non-Vented
without Filter
700032



Drip Chamber Inline Non-Vented
with 15 Micron Filter
700547001



Drip Chamber Non-Filtered
700442001

MOLDED PARTS & COMPONENTS

Drip Chambers

Cap: Polyethylene Spike/Inline Tubing Adapter: ABS
Chamber: Polypropylene Vent: Polyethylene



Drip Chamber, Vented without Filter

700442002



Drip Chamber, Vented without Filter, Large Bore Tube Connector

701242001



Drip Chamber Non-Vented Bent Stem

700541001



Drip Chamber Non-Vented Filter

700029



Drip Vented with Ball Valve
Ball: Polypropylene
Chamber: PVC

701610001

MOLDED PARTS &
COMPONENTS