


## Fixation Devices


The Revolution™ is a catheter securement device that holds nonvascular catheters in place for up to seven days without sutures. The Revolution helps minimize catheter movement and accidental removal. It includes a lid to keep the site dry when showering.



### FEATURES & BENEFITS


Wide range	Secures nonvascular catheters from 6F-24F
Central adhesive hydrocolloid	Skin-friendly, absorbent around puncture site, and conforms well to body anatomy
Advanced adhesive	Provides safe and secure fixation for up to seven days
Included showercap	Covers insertion site for added protection while showering

## Fixation Devices



Revolution™ Needleless  
Non Vascular Catheter Securement  
Device (6F-24F)

**CFG624**



StayFIX® Securement Device for  
Percutaneous Catheters

5F-14F **680ME**

14F-22F **685ME**