

## Roller Clamps



Small Bore Roller Clamp

**700022001**



Large Bore Roller Clamp

**700022002**

## Drip Chambers

Many tubing sets available with non-DEHP material.

Cap: Polyethylene Spike/Inline Tubing Adapter: ABS Chamber: Polypropylene Vent: Polyethylene



Macro Inline Non-Vented

**700098**



Drip Chamber, Vented with 15 Micron Filter

**700026**



Drip Chamber Non-Vented without Filter

**700032**



Drip Chamber Inline Non-Vented with 15 Micron Filter

**700547001**