













# Extension Tubing

High Durometer, Non Expanding.

Material: PVC, Clear non-DEHP

						
<b>Male to Male</b>	6" Length (15cm)	12" Length (30cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)
	<b>FL6006P</b>	<b>FL6012P</b>	<b>FL6024P</b>	<b>FL6036P</b>	<b>FL6048P</b>	<b>FL6060P</b>
						
<b>Male to Female</b>	6" Length (15cm)	12" Length (30cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)
	<b>FL6106P</b>	<b>FL6112P</b>	<b>FL6124P</b>	<b>FL6136P</b>	<b>FL6148P</b>	<b>FL6148CP*</b> <b>FL6160P</b>

\*unique coil

# Fluid Administration Roller Clamp



Small Bore Roller Clamp

**700022001**



Large Bore Roller Clamp

**700022002**