

Fluid Administration Sets

Small Bore with Drip Chamber. Many tubing sets available with non-DEHP material

Tubing Material: PVC



48" Small Bore, Vented Macro Drip Chamber without Micro Filter, Male Luer ID: .120", Inches: 48", OD: .160"

FAS1048P



60" Small Bore, Vented Macro Drip Chamber, without Micro Filter, Male Luer ID: .120", Inches: 60", OD: .160"

FAS1060P



72" Small Bore, Vented Macro Drip Chamber, without Micro Filter, Male Luer ID: .120", Inches: 72", OD: .160"

FAS1072P



108" Small Bore, Vented Macro Drip Chamber, with Micro Filter, Male Luer ID: .120", Inches: 60", OD: .160"

FAS1108P

Fluid Administration Sets

Small Bore with Drip Chamber

Tubing Material: PVC



60" Small Bore, Vented Macro Drip Chamber without Micro Filter, Male Luer
ID: .120", Inches: 60", OD: .160"

FAS1160P



72" Small Bore, Vented Macro Drip Chamber without Micro Filter, Male Luer
ID: .120", Inches: 72", OD: .160"

FAS1172P



60" Small Bore, Non-Vented Macro Drip Chamber without Micro Filter, Male Luer
ID: .120", Inches: 60", OD: .160"

FAS1260P



72" Small Bore, Non-Vented Macro Drip Chamber without Micro Filter, Male Luer
ID: .120", Inches: 72", OD: .160"

FAS1272P

FLUID
MANAGEMENT

Fluid Administration Sets

Small Bore with Drip Chamber

Tubing Material: PVC



24" Small Bore with User Selectable Macro Drip Chamber without Micro Filter, 48" Small Bore Bifurcated with Roller Clamps, Male Luers ID: .120", Inches: 24", OD: .160"

FAS2072P



24" Small Bore with User Selectable Macro Drip Chamber without Micro Filter, 48" Small Bore Trifurcated with Roller Clamps, Male Luers ID: .120", Inches: 24", OD: .160"

FAS3072P



48" Small Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer ID: .120", Inches: 60", OD: .160"

FAS4048P




60" Small Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer ID: .120", Inches: 60", OD: .160"

FAS4060P

Fluid Administration Sets

Small and Large Bore with Drip Chamber

Tubing Material: PVC



72" Small Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer
ID: .120", Inches: 72", OD: .160"

FAS4072P



108" Large Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer
ID: .120", Inches: 108", OD: .160"

FAS4108P



48" Large Bore with User Selectable Vented Macro Drip Chamber with Micro Filter, 48" Large Bore, Bifurcated with Roller Clamps, Male Luers
ID: .120", Inches: 48", OD: .160"

FAS5048P



50" Large Bore with User Selectable Vented Macro Drip Chamber with Micro Filter, 50" Large Bore, Bifurcated with Roller Clamps, Male Luers
ID: .120", Inches: 50", OD: .160"

FAS5060P

Fluid Administration Sets

Small and Large Bore with Drip Chamber

Tubing Material: PVC



72" Large Bore with User Selectable Vented Macro Drip Chamber with Micro Filter, 72" Large Bore Bifurcated with Roller Clamps, Male Luers
ID: .72", Inches: 50", OD: .160"

FAS5072P



108" Large Bore with User Selectable Vented Macro Drip Chamber with Micro Filter, 108" Large Bore, Bifurcated with Roller Clamps, Male Luers
ID: .108", Inches: 108", OD: .160"

FAS5108P



48" Large Bore Non-vented Macro Drip Chamber with Micro Filter, 48" Large Bore, with Roller Clamps, Male Luers
ID: .120", Inches: 48", OD: .160"

FAS5248P



24" Small Bore with Vented Macro Drip Chamber with Micro Filter, 24" Small Bore Bifurcated with Roller Clamps, Male Luers
ID: .24", Inches: 524", OD: .160"

FAS6072

FLUID
MANAGEMENT

Fluid Administration Sets

Large Bore with Drip Chamber

Tubing Material: PVC



72" Large Bore, Piercing Spike with Sosa Vented Macro Drip Chamber with Micro Filter, Male Luer ID: .120", Inches: 72", OD: .160"

K08-00006



2" Large Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer ID: .120", Inches: 724", OD: .160", 3-Way Stopcock

K08-00077B



2" Large Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer ID: .120", Inches: 724", OD: .160", 1-Way Stopcock

K08-00078A



6.5" Large Bore, Vented Macro Drip Chamber with Micro Filter, Male Luer, 3-Way Stopcock ID: .120", Inches: 72", OD: .160"

K08-00417BP



6.5" Large Bore, Vented Spike Macro Drip Chamber with Micro Filter, Male Luer, 3-Way Stopcock ID: .120", Inches: 72", OD: .160"

K08-00577

FLUID
MANAGEMENT

Fluid Administration Sets

Large Bore with Drip Chamber

Tubing Material: PVC



6.5" Large Bore, Piercing Spike with Sosa Vented Macro Drip Chamber with Micro Filter, Male Luer
ID: .120", Inches: 72", OD: .160" 3-way Stopcock

K08-00685A



6.5" Large Bore, Vented Spike Macro Drip Chamber with Micro Filter, Male Luer, 3-Way Stopcock with Swivel Nut
ID: .6.4", Inches: 72", OD: .160"

K08-01878B