



ConcierGE® Guiding Catheters

Our advanced guiding catheters are engineered with braiding technology, distal segments, and a proprietary polymer-blend shaft for exceptional back up support.



1:1 TORQUE

The ConcierGE demonstrates excellent 1:1 torque response, promoting quicker catheter seating with less whipping.



BACKUP SUPPORT

The ConcierGE provides reliable backup support with proprietary nylon blend shaft material and wire braid technology giving you the ability to advance balloons and stents to the diseased artery without losing shaft strength.

ATRAUMATIC TIP

The ConcierGE has one of the softest, most radiused tips for safer tracking and seating so physicians of all experience levels can navigate challenging anatomy with confidence.



KINK RESISTANCE


The ConcierGE has approximately 30% better kink resistance* than the leading brands, facilitating the advancement and withdrawal of interventional devices.

* Test Results are based on 6F guide catheter comparison. Data on file.




JUDKINS RIGHT

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100JR30	5	3.0	-	CGC7100JR30	7	3.0	-
CGC5100JR35	5	3.5	-	CGC7100JR35	7	3.5	-
CGC5100JR35SH	5	3.5	2	CGC7100JR35SH	7	3.5	2
CGC5100JR40	5	4.0	-	CGC7100JR40	7	4.0	-
CGC5100JR40SH	5	4.0	2	CGC7100JR40SH	7	4.0	2
CGC5100JR50	5	5.0	-	CGC7100JR50	7	5.0	-
CGC5100JR50SH	5	5.0	2	CGC7100JR50SH	7	5.0	2
CGC5100JR60	5	6.0	-	CGC7100JR60	7	6.0	-
CGC6100JR30	6	3.0	-	CGC8100JR30	8	3.0	-
CGC6100JR35	6	3.5	-	CGC8100JR35	8	3.5	-
CGC6100JR35SH	6	3.5	2	CGC8100JR35SH	8	3.5	2
CGC6100JR40	6	4.0	-	CGC8100JR40	8	4.0	-
CGC6100JR40SH	6	4.0	2	CGC8100JR40SH	8	4.0	2
CGC6100JR50	6	5.0	-	CGC8100JR50	8	5.0	-
CGC6100JR50SH	6	5.0	2	CGC8100JR50SH	8	5.0	2
CGC6100JR60	6	6.0	-	CGC8100JR60	8	6.0	-



JUDKINS LEFT

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100JL30	5	3.0	-	CGC7100JL30	7	3.0	-
CGC5100JL35	5	3.5	-	CGC7100JL35	7	3.5	-
CGC5100JL35SH	5	3.5	2	CGC7100JL35SH	7	3.5	2
CGC5100JL40	5	4.0	-	CGC7100JL40	7	4.0	-
CGC5100JL40SH	5	4.0	2	CGC7100JL40SH	7	4.0	2
CGC5100JL45	5	4.5	-	CGC7100JL45	7	4.5	-
CGC5100JL45SH	5	4.5	2	CGC7100JL45SH	7	4.5	2
CGC5100JL50	5	5.0	-	CGC7100JL50	7	5.0	-
CGC5100JL50SH	5	5.0	2	CGC7100JL50SH	7	5.0	2
CGC6100JL30	6	3.0	-	CGC8100JL30	8	3.0	-
CGC6100JL35	6	3.5	-	CGC8100JL35	8	3.5	-
CGC6100JL35SH	6	3.5	2	CGC8100JL35SH	8	3.5	2
CGC6100JL40	6	4.0	-	CGC8100JL40	8	4.0	-
CGC6100JL40SH	6	4.0	2	CGC8100JL40SH	8	4.0	2
CGC6100JL45	6	4.5	-	CGC8100JL45	8	4.5	-
CGC6100JL45SH	6	4.5	2	CGC8100JL45SH	8	4.5	2
CGC6100JL50	6	5.0	-	CGC8100JL50	8	5.0	-
CGC6100JL50SH	6	5.0	2	CGC8100JL50SH	8	5.0	2



ConcierGE® Guiding Catheters

JUDKINS LEFT ST

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100JL35ST	5	3.5	-	CGC7100JL35ST	7	3.5	-
CGC5100JL40ST	5	4.0	-	CGC7100JL40ST	7	4.0	-
CGC6100JL35ST	6	3.5	-	CGC8100JL35ST	8	3.5	-
CGC6100JL40ST	6	4.0	-	CGC8100JL40ST	8	4.0	-



SBS30

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100SBS30	5	3.0	-	CGC7100SBS30	7	3.0	-
CGC5100SBS30SH	5	3.0	2	CGC7100SBS30SH	7	3.0	2
CGC6100SBS30	6	3.0	-	CGC8100SBS30	8	3.0	-
CGC6100SBS30SH	6	3.0	2	CGC8100SBS30SH	6	3.0	2



SBS35

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100SBS35	5	3.5	-	CGC7100SBS35	7	3.5	-
CGC5100SBS35SH	5	3.5	2	CGC7100SBS35SH	7	3.5	2
CGC6100SBS35	6	3.5	-	CGC8100SBS35	8	3.5	-
CGC6100SBS35SH	6	3.5	2	CGC8100SBS35SH	8	3.5	2



SBS 4

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100SBS40	5	4.0	-	CGC7100SBS40	7	4.0	-
CGC5100SBS40SH	5	4.0	2	CGC7100SBS40SH	7	4.0	2
CGC6100SBS40	6	4.0	-	CGC8100SBS40	8	4.0	-
CGC6100SBS40SH	6	4.0	2	CGC8100SBS40SH	8	4.0	2



ConcierGE® Guiding Catheters

AMPLATZ RIGHT 1

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AR1	5	AR1	-	CGC7100AR1	7	AR1	-
CGC5100AR1SH	5	AR1	2	CGC7100AR1SH	7	AR1	2
CGC6100AR1	6	AR1	-	CGC8100AR1	8	AR1	-
CGC6100AR1SH	6	AR1	2	CGC8100AR1SH	8	AR1	2



AMPLATZ RIGHT 2

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AR2	5	AR2	-	CGC7100AR2	7	AR2	-
CGC5100AR2SH	5	AR2	2	CGC7100AR2SH	7	AR2	2
CGC6100AR2	6	AR2	-	CGC8100AR2	8	AR2	-
CGC6100AR2SH	6	AR2	2	CGC8100AR2SH	8	AR2	2




AMPLATZ RIGHT 3

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AR3	5	AR3	-	CGC7100AR3	7	AR3	-
CGC6100AR3	6	AR3	-	CGC8100AR3	8	AR3	-



AMPLATZ LEFT 1

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AL1	5	AL1	-	CGC7100AL1	7	AL1	-
CGC5100AL1SH	5	AL1	2	CGC7100AL1SH	7	AL1	2
CGC6100AL1	6	AL1	-	CGC8100AL1	8	AL1	-
CGC6100AL1SH	6	AL1	2	CGC8100AL1SH	8	AL1	2



ConcierGE® Guiding Catheters

AMPLATZ LEFT 2

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AL2	5	AL2	-	CGC7100AL2	7	AL2	-
CGC5100AL2SH	5	AL2	2	CGC7100AL2SH	7	AL2	2
CGC6100AL2	6	AL2	-	CGC8100AL2	8	AL2	-
CGC6100AL2SH	6	AL2	2	CGC8100AL2SH	8	AL2	2

AMPLATZ LEFT 2

AMPLATZ LEFT 3

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100AL3	5	AL3	-	CGC7100AL3	7	AL3	-
CGC6100AL3	6	AL3	-	CGC8100AL3	8	AL3	-

AMPLATZ LEFT 3

MULTIPURPOSE 1

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100MP1	5	MP	-	CGC7100MP1	7	MP	-
CGC5100MP1SH	5	MP	2	CGC7100MP1SH	7	MP	2
CGC6100MP1	6	MP	-	CGC8100MP1	8	MP	-
CGC6100MP1SH	6	MP	2	CGC8100MP1SH	8	MP	2

MULTIPURPOSE 1

LEFT CORONARY BYPASS

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100LCB	5	LCB	-	CGC7100LCB	7	LCB	-
CGC5100LCBSH	5	LCB	2	CGC5100LCBSH	5	LCB	2
CGC6100LCB	6	LCB	-	CGC8100LCB	8	LCB	-
CGC6100LCBSH	6	LCB	2	CGC8100LCBSH	8	LCB	2

LEFT CORONARY BYPASS

CATHETERS & EXTRUSIONS

ConcierGE® Guiding Catheters

RIGHT CORONARY BYPASS

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100RCB	5	RCB	-	CGC7100RCB	7	RCB	-
CGC5100RCBSH	5	RCB	2	CGC5100RCBSH	5	RCB	2
CGC6100RCB	6	RCB	-	CGC8100RCB	8	RCB	-
CGC6100RCBSH	6	RCB	2	CGC8100RCBSH	8	RCB	2



INTERNAL MAMMARY

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100IM	5	IM	-	CGC7100IM	7	IM	-
CGC5100IMSH	5	IM	2	CGC7100IMSH	7	IM	2
CGC6100IM	6	IM	-	CGC8100IM	8	IM	-
CGC6100IMSH	6	IM	2	CGC8100IMSH	8	IM	2



EL GAMAL

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100EG	5	EG	-	CGC7100EG	7	EG	-
CGC6100EG	6	EG	-	CGC8100EG	8	EG	-



HOCKEY STICK

Catalog Number	French Size	Tip	Side Holes	Catalog Number	French Size	Tip	Side Holes
CGC5100HS	5	HS	-	CGC7100HS	7	HS	-
CGC6100HS	6	HS	-	CGC8100HS	8	HS	-

