

Hemostasis Valves

AccessPLUS™

Material: Polycarbonate Inner Lumen: 9F (.120")
Valve Material: Silicone Luer: Rotating



Large Bore Hemostasis Valve
500628002



Large Bore Hemostasis Valve with
12" PVC Tubing Extension
563010001



Large Bore Hemostasis Valve with
6" Polyurethane Tubing Extension
563009001

DoublePlay™

Material: Polycarbonate Inner Lumen: 9F (.120")
Valve Material: Silicone Luer: Rotating



French Size: 9F, Tubing Extension:
No Tubing Extension, Tubing Material: No
Tubing Material, Brand: DoublePlay™

501075001

Hemostasis Valves

Honor®

Material: Polycarbonate Inner Lumen: 7F (.094")
Valve Material: Silicone Luer: Rotating



Hemostasis Valve
580195001



Hemostasis Valve with
8" Tubing Extension
MAP350

MBA™

Material: Polycarbonate Inner Lumen: 9F (.118")
Valve Material: Silicone Luer: Rotating



Large Bore Dual Seal
Hemostasis Valve
500981002



Hemostasis Valve with 10"
PVC Tubing Extension
501086001

Hemostasis Valves

Passage®

Material: Polycarbonate
Valve Material: Silicone



Hemostasis Valve with
Rotating Luer

500066



Hemostasis Valve with Slip Fit
(.179" to .209" taper)
Inner Lumen: 7.5F (.097")

K12-00457P