

Custom Tubing

Over the last 30 years, Merit has developed the infrastructure needed for mass customization to meet the unique needs of individual clinics or physicians. This translates into a robust OEM service that allows you to receive a customized configuration as a bulk non-sterile product or even a private labeled sterile product.



TUBING

ASPIRATION SETS • TUBING EXTENSION SETS • FLUID ADMINISTRATION SETS
HIGH PRESSURE TUBING • PRESSURE MONITORING TUBING



HOW MERIT HELPS OUR OEM CUSTOMERS:

- Tubing Set Design and Development
- Inline Quality Control
- Packaging Validations
- Sterilization Validations
- Private Labeling
- Regulatory Support






High Pressure Tubing

Rated at 1200 PSI (83 BAR)

Material: Flexible Nylon Braided Polyurethane

Male to Female

Airless Rotating Adapter

				
10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	72" Length (183cm)
HPF100P	HPF200P	HPF300P	HPF480P	HPF720P

Fixed Male Connector

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	72" Length (183cm)
HPF101P	HPF201P	HPF301P	HPF481P	HPF721P

High Pressure Tubing






Rated at 1200, 900 and 500 PSI (83, 62 and 34 BAR)

Material: PVC, Clear

Male to Female

Airless Rotating Adapter

1200 PSI (83 BAR)
900 PSI (62 BAR)
500 PSI (34 BAR)

				
10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	59" Length (150cm)
HPT100P HP9100P HP5100P	HPT200P HP9200P HP5200P	HPT300P HP9300P HP5300P	HPT480P HP9480P HP5480P	HPT590P

Fixed Male Connector

1200 PSI (83 BAR)
900 PSI (62 BAR)
500 PSI (34 BAR)

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	48" Length (122cm)	59" Length (150cm)
HPT101P HP9101P HP5101P	HPT201P HP9201P HP5201P	HPT301P HP9301P HP5301P	HPT481P HP9481P HP5481P	HPT591P




Excite™ High Pressure Tubing

Rated at 1200 PSI (83 BAR)

Material: Coextruded Nylon/Polyurethane

TUBING

Male to Female



Airless Rotating Adapter

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	39" Length (100cm)	48" Length (122cm)	59" Length (150cm)
HPC100P	HPC200P	HPC300P	HPC390P	HPC480P	HPC590P

Fixed Male Connector

10" Length (25cm)	20" Length (51cm)	30" Length (76cm)	39" Length (100cm)	48" Length (122cm)	59" Length (150cm)
HPC101P	HPC201P	HPC301P	HPC391P	HPC481P	HPC591P

