

Catheter Hubs

Material: Polycarbonate



4F Hub, White
100941200



5F Hub, White
100941205



6F Hub, White
.050" thru-hole
100941211




6F Hub, White
100941212



8F Hub, White
100941218

Hub details. Measurements in inches.



C $^{+.002}_{-.001}$ \varnothing at entrance straightline tapers to **D** $^{+.001}_{-.001}$ \varnothing at bottom

Catalog Number	A	B	C	D
100941200	0.087	0.090	0.058	0.056
100941205	0.111	0.090	0.071	0.069
100941211	0.120	0.125	0.083	0.081
100941212	0.012	0.0125	0.083	0.081
100941218	0.135	0.125	0.109	0.081



Custom colors
and printing
options available.



Reverse

5F Hub, White
0.045 inner lumen
Printed

590012039



Reverse

6F Hub, White
0.050 inner lumen
Printed

590012057