



Inflation Devices

- A worldwide market leader
- Ergonomic handle design
- Pressures up to 35ATM (514PSI)
- 5ml of safety space intended to minimize air bubbles introduced into the balloon
- Packaged stand alone or customized to your specifications
- Digital technologies

SEE RELATED PRODUCTS



VacLok® Syringes
220



Stopcocks
253



Hemostasis Valves
242

May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.

Blue Diamond™ Digital Inflation Device

Combines the accuracy of a digital device with the simplicity of an analog system. Features one button functionality for ease of use and a large dial face with oversize time and pressure digits to enhance visibility of critical information. Available stand alone or call to customize to your specifications.

Display: Digital

Pressure: 30 ATM

Stopcock: 3-Way 500 PSI (34 Bar)

Tubing: 13" (33.02cm) Extension Tubing with Male Rotating Luer

Barrel

Volume: 20ml



- Timed Inflation Duration
- Timed Elapsed Since Last Inflation
- Displays pressures in ATM/BAR or PSI

Blue Diamond™ Inflation Device

IN7130

BasixCOMPAK™ Inflation Device

Ergonomically designed to ensure accurate inflation of interventional devices. Features an angled analog gauge for better viewing of pressure readings and a luminescent dial for low-lit procedures. Available stand alone or call to customize to your specifications.

Display: Analog

Pressure: 30 ATM

Stopcock: 3-Way 500 PSI (34 Bar)

Tubing: 13" (33.02cm) Extension Tubing with Male Rotating Luer

Barrel

Volume: 20ml



BasixCOMPAK™ Inflation Device

IN4130

Limited™ Inflation Device

No pressure display. Available stand alone or call to customize to your specifications.

Pressure: -1 to 25 ATM/BAR, 0 to 367 PSI

Tubing: 13" (33.02cm) Extension Tubing with Male Rotating Luer

Barrel

Volume: 20ml



Limited™ Inflation Device

K05-01299

Big60™ Inflation Device

Ergonomically designed to ensure accurate inflation and deflation of interventional devices. Features an angled analog gauge for better viewing of pressure readings and a luminescent dial for low-lit procedures. Designed to monitor inflation pressures over a range of -.68 to 12 ATM/BAR (-10 to 176 psi) with an accuracy of 3% ATM/BAR. Available stand alone or call to customize to your specifications.

Display: Analog

Pressure: 12 ATM/BAR

Stopcock: 3-Way

Tubing: 13" (33.02cm) Extension Tubing with Male Rotating Luer

Barrel Volume: 60ml



Big60™ Inflation Device

ENDO-AN6012

basixTOUCH™

The basixTOUCH delivers the speed you need to expand the capacity of your lab. Quick-release handle for easy release even from high-pressure inflations, accelerating deflation time. Ergonomic handle and new threads ensure that high pressures can be reached rapidly – up to 35 atm with less than 3 rotations of the handle. One-handed preparation.

Display: Analog

Pressure: 35 ATM/BAR

Tubing: 13" (33.02cm) Extension Tubing with Male Rotating Luer

Barrel
Volume: 30ml



basixTOUCH™ Inflation Device

IN8130