

# Coaxial Introducer/Dilator Pairs



French Size: 4F  
Length: 10cm

**MAK4F-10**

French Size: 5F  
Length: 10cm

**MAK5F-10**

French Size: 5F  
Length: 15cm

**MAK5F-15**

French Size: 4F  
Length: 10cm

**S-MAK4F-10**  
Stiffened Dilator

French Size: 5F  
Length: 10cm

**S-MAK5F-10**  
Stiffened Dilator

French Size: 5F  
Length: 15cm

**S-MAK5F-15**  
Stiffened Dilator

PERCUTANEOUS  
ACCESS

# Coaxial Introducers



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK401**



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK401-E**



French Size: 4F  
 Needle Length: 4cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK401-PED**



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm  
 Gauge: 21, Guide Wire Diameter/Length:  
 0.018" (0.46mm) x 40cm

**MAK401N**



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm  
 Gauge: 21, Guide Wire Diameter/Length:  
 0.018" (0.46mm) x 40cm, Echo Enhanced

**MAK401N-E**



French Size: 4F  
 Needle Length: 4cm  
 Core/Tip: Nitinol/Platinum Tip, Length: 10cm  
 Gauge: 21, Guide Wire Diameter/Length:  
 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK401N-PED**

PERCUTANEOUS  
ACCESS

# Coaxial Introducers



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Platinum Tip  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm

**MAK405**



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Platinum Tip  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK405-E**



French Size: 4F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK405-PED**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21,  
 Guide Wire Diameter/Length: 0.018"  
 (0.46mm) x 40cm

**MAK501**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK501-E**



French Size: 5F  
 Needle Length: 4cm  
 Core/Tip: Stainless Steel/Stainless Steel  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK501-PED**

# Coaxial Introducers



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Nitinol/Platinum Tip  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm

**MAK501N**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Nitinol/Platinum Tip  
 Length: 15cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm

**MAK501N-15**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Nitinol/Platinum Tip  
 Length: 10cm, Gauge: 21, Guide Wire  
 Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK501N-E**

PERCUTANEOUS ACCESS



French Size: 5F  
 Needle Length: 4cm  
 Core/Tip: Stainless Steel/Stainless Steel, Length: 10cm, Gauge: 21  
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK501N-PED**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21  
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm

**MAK505**



French Size: 5F  
 Needle Length: 7cm  
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21  
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK505-E**



French Size: 5F  
 Needle Length: 4cm  
 Core/Tip: Stainless Steel/Platinum Tip, Length: 10cm, Gauge: 21  
 Guide Wire Diameter/Length: 0.018" (0.46mm) x 40cm  
 Echo Enhanced

**MAK505-PED**

## Guide Sheaths

Contact your Merit Medical OEM representative for more information about our guide sheath options.



PERCUTANEOUS  
ACCESS

# Large Bore Introducers



French Size: 16F  
Sheath Length: 13cm  
Dilator Length: 20cm  
With Sideport

**FCL-106-00**



French Size: 16F  
Sheath Length: 13cm  
Dilator Length: 20cm  
No Sideport

**FCL-108-00**

PERCUTANEOUS  
ACCESS

## Long Introducers

Merit's Long Introducers have a braided core for excellent torqueability. The tapered dilator and soft, atraumatic tip is designed for enhanced tissue crossing. The embedded platinum marker aids in radiopacity for proper placement.



Catalog Number	French Size	Length	Curve Angle
<b>FCL-160-00</b>	8.5F	60cm	15° Curve
<b>FCL-160-01</b>	8.5F	60cm	30° Curve
<b>FCL-160-02</b>	8.5F	60cm	55° Curve
<b>FCL-160-03</b>	8.5F	60cm	120° Curve (short)
<b>FCL-160-04</b>	8.5F	60cm	120° Curve (long)
<b>FCL-160-05</b>	8.5F	60cm	150° Curve
<b>FCL-160-06</b>	8.5F	79.4cm	15° Curve
<b>FCL-160-07</b>	8.5F	79.4cm	55° Curve
<b>FCL-160-08</b>	8.5F	79.4cm	90° Curve
<b>FCL-160-09</b>	8.5F	79.4cm	120° Curve
<b>FCL-160-10</b>	8.5F	101.5cm	55° Curve