

Long Introducers

Merit's Long Introducers have a braided core for excellent torqueability. The tapered dilator and soft, atraumatic tip is designed for enhanced tissue crossing. The embedded platinum marker aids in radiopacity for proper placement.



Catalog Number	French Size	Length	Curve Angle
FCL-160-00	8.5F	60cm	15° Curve
FCL-160-01	8.5F	60cm	30° Curve
FCL-160-02	8.5F	60cm	55° Curve
FCL-160-03	8.5F	60cm	120° Curve (short)
FCL-160-04	8.5F	60cm	120° Curve (long)
FCL-160-05	8.5F	60cm	150° Curve
FCL-160-06	8.5F	79.4cm	15° Curve
FCL-160-07	8.5F	79.4cm	55° Curve
FCL-160-08	8.5F	79.4cm	90° Curve
FCL-160-09	8.5F	79.4cm	120° Curve
FCL-160-10	8.5F	101.5cm	55° Curve