

Luers & Connectors

SEE RELATED
PRODUCTS



Caps & Covers
103



Tubing
223

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

Female Slip Fit Luers

Material: Polycarbonate



White Female Luer Lock
Slip Fit .120" to .150" Taper
100916001



Clear Female Luer Lock
Slip Fit .120" to .150" Taper
102112001



Clear Female Luer Lock
Slip Fit .154" to .192" Taper
102102001



Gamma Color Female Luer Lock
Slip Fit .154" to .192" Taper
102202001



Gamma Color Female Luer Lock
Slip Fit .120" to .150" Taper
102202002



Clear Female Luer Lock
Slip Fit .179" to .209" Taper
100719

Female Tubing Pocket Fit Luers

Material: Polycarbonate



Clear Female Luer Lock
Tubing Pocket .095" to .101"

102103002



Royal Blue Female Luer Lock
Tubing Pocket .095" to .101"

100972002



Clear Female Luer Lock
Tubing Pocket .118" to .124"

101023001



Yellow Female Luer Lock
Tubing Pocket .118" to .124"

101023002



Clear Female Luer Lock
Tubing Pocket .124" to .132"

102103001



Clear Female Luer Lock
Tubing Pocket .140" to .145"

102103004



Clear Female Luer Lock
Tubing Pocket .158" to .164"

102103003



Clear Female Luer Lock
Tubing Pocket .182" to .192"

100976001



Gamma Color Female Luer Lock
Tubing Pocket .095" to .101"

102203002



Gamma Color Female Luer Lock
Tubing Pocket .118" to .124"

101123001



Gamma Color Female Luer Lock
Tubing Pocket .124" to .132"

102203001



Gamma Color Female Luer Lock
Tubing Pocket .158" to .164"

102203003



Clear Female Luer Lock
Tubing Pocket .124" to .132"

101071002

MOLDED PARTS & COMPONENTS

Male Slip Fit Luers

Material: Polycarbonate



Blue Male Luer Lock
Slip Fit .120" to .150" Taper
100914002



Gamma Color Male Luer Lock
Slip Fit .120" to .150" Taper
102200002



Gamma Color Male Luer Lock
Slip Fit .154" to .192" Taper
102200001



White Male Luer Lock
Slip Fit .120" to .150" Taper
100914001



Clear Male Luer Lock
Slip Fit .120" to .150" Taper
102110001



Clear Male Luer Lock
Slip Fit .154" to .192"
102100001



Gamma Color, Male Luer
Lock Slip Fit .148" to .166"
102056002



Gamma Color Male Luer Lock
Slip Fit .178" to .208"
102039002



Clear Male Luer Lock
Slip Fit .148" to .166"
102056001



Clear Male Luer Lock
Slip Fit .178" to .208"
102039

MOLDED PARTS &
COMPONENTS

Male Tubing Pocket Fit Luers

Material: Polycarbonate



Clear Male Luer Lock
Tubing Pocket .095" to .101"

102101002



Clear Male Luer Lock
Tubing Pocket .108" to .114"

102101005



Clear Male Luer Lock
Tubing Pocket .118" to .124"

101024001



Clear Male Luer Lock
Tubing Pocket .124" to .132"

102101001



Clear Male Luer Lock
Tubing Pocket .140" to .145"

102101004



Clear Male Luer Lock
Tubing Pocket .158" to .164"

102101003



Clear Male Luer Lock
Tubing Pocket .182" to .192"

100975001



Gamma Color Male Luer Lock
Tubing Pocket .118" to .124"

101124001



Gamma Color Male Luer Lock
Tubing Pocket .124" to .132"

102201001



Gamma Color Male Luer Lock
Tubing Pocket .095" to .101"

102201002



Gamma Color Male Luer Lock
Tubing Pocket .158" to .164"

102201003

MOLDED PARTS &
COMPONENTS

Luers & Connectors Other



Clear, 2-Way Connector
Tubing Pockets
Connector 1: .160" Tubing Pocket
Connector 2: .160" Tubing Pocket
Connector 3: .160" Tubing Pocket
Material: Cyrolite

101002



Clear, 3-Way Connector
Tubing Pockets
Connector 1: .160" Tubing Pocket
Connector 2: .160" Tubing Pocket
Connector 3: .160" Tubing Pocket
Material: Cyrolite

101003



Y-Connector
Material: PVC

700155002



Clear Male/Female Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Female Luer Lock
Material: Polycarbonate

101047001



Red Male/Female Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Female Luer Lock
Material: Polycarbonate

101542001



Blue Male/Female Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Female Luer Lock
Material: Polycarbonate

101542002



Clear Female/Female Luer
Lock Connector
Connector 1: Female Luer Lock
Connector 2: Female Luer Lock
Material: MABS

700357001



Clear, 2-Way Connector
Tubing Pockets
Connector 1: .160" Tubing Pocket
Connector 2: .160" Tubing Pocket
Connector 3: .160" Tubing Pocket
Material: Polycarbonate

101004



Y Tubing Connector
Connector 1: .154 slip
Connector 2: .154 slip
Connector 3: .154 slip
Material: Clear ABS

101508001



Clear Male/Male Luer
Lock Connector
Connector 1: Male Luer Lock
Connector 2: Male Luer Lock
Material: Polycarbonate

101025001

T-Connectors

Material: Polycarbonate



Clear Two Female Luer Lock,
One Male Slip

100966001



Gamma Color Two Female
Luer Lock, One Male Slip

100966002



Clear One Female Luer Lock,
Two Male Slip

100969001



Clear All Male Slip

100967001



Clear Two Female Luer Lock,
One Male Swivel

501314001



Clear All Female Luer Lock

100965001



Clear with .139" Tubing Pockets

100965002