




Microcatheters

- Flexible (distal 200mm) design for atraumatic vessel entry.
- Hydrophilic coating (distal 80cm) for smooth traction through tortuous anatomy.
- New swan neck design and high performance of preshaped distal end are ideal for access of hepatic artery.
- Radiopaque marker for easy identification of distal tip of catheter. Located 1.3mm proximal to the distal tip.



Why Use Merit Maestro® Microcatheters?

- All Merit Maestro® Microcatheters are rated to 600psi
- Embolics up to 700µm (microns)
- Coils .018"

RELATED PRODUCTS



Guide Wires



Needles

May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.

Merit Maestro® Microcatheters 2.8F Tapered to 2.4F

Catalog Number	Length		Length (cm)	Tip	Minimum Inner Diameter	Guide Wire	Guide Catheter Indicated for Guide Wire Size
	Proximal	Distal					
28MC24110ST	2.8F	2.4F	110	Straight	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC24130ST	2.8F	2.4F	130	Straight	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC24150ST	2.8F	2.4F	150	Straight	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC2411045	2.8F	2.4F	110	45 degree	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC2413045	2.8F	2.4F	130	45 degree	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC2415045	2.8F	2.4F	150	45 degree	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC24110SN	2.8F	2.4F	110	Swan Neck	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC24130SN	2.8F	2.4F	130	Swan Neck	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)
28MC24150SN	2.8F	2.4F	150	Swan Neck	.020" (0.53mm)	.018" (0.46mm)	.038" (0.97mm)

Merit Maestro® Microcatheters 2.8 No Taper

Catalog Number	Length		Length (cm)	Tip	Minimum Inner Diameter	Guide Wire	Guide Catheter Indicated for Guide Wire Size
	Proximal	Distal					
28MC28110ST	2.8F	2.8F	110	Straight	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC28130ST	2.8F	2.8F	130	Straight	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC28150ST	2.8F	2.8F	150	Straight	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC2811045	2.8F	2.8F	110	45 degree	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC2813045	2.8F	2.8F	130	45 degree	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC2815045	2.8F	2.8F	150	45 degree	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC28110SN	2.8F	2.8F	110	Swan Neck	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC28130SN	2.8F	2.8F	130	Swan Neck	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)
28MC28150SN	2.8F	2.8F	150	Swan Neck	.024" (0.64mm)	.021" (0.53mm)	.038" (0.97mm)

Merit Maestro® Microcatheters 2.9 No Taper

Catalog Number	Length		Length (cm)	Tip	Minimum Inner Diameter	Guide Wire	Guide Catheter Indicated for Guide Wire Size
	Proximal	Distal					
29MC29110ST	2.9F	2.9F	110	Straight	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC29130ST	2.9F	2.9F	130	Straight	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC29150ST	2.9F	2.9F	150	Straight	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC2911045	2.9F	2.9F	110	45 degree	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC2913045	2.9F	2.9F	130	45 degree	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC2915045	2.9F	2.9F	150	45 degree	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC29110SN	2.9F	2.9F	110	Swan Neck	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC29130SN	2.9F	2.9F	130	Swan Neck	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)
29MC29150SN	2.9F	2.9F	150	Swan Neck	.027" (0.68mm)	.021" (0.53mm)	.042" (1.07mm)