



# Sheath Introducers

- 135 Prelude® Sheath Introducers 11cm
- 136 Prelude® Sheath Introducers 23cm
- 137 Prelude PRO™ Sheath Introducers
- 137 Prelude ACT™ Sheath Introducers
- 138 Prelude® Sheath Introducers with Merit Advance® Needle
- 139 Prelude® Sheath Introducers with Marker Tip
- 140 Prelude® Sheath Introducers – Mini Access Kits
- 141 Prelude® Short Sheath Introducers

SEE RELATED  
PRODUCTS



Needles  
145




Guide Wires  
89



Catheters  
5

May not be available for all OEM applications.

# Prelude® Sheath Introducers & Dilators



PERCUTANEOUS ACCESS

# Prelude® Sheath Introducers 11cm

The Prelude® Sheath Introducer is packaged in a holster that keeps all components secure, especially the guide wire. Components can be flushed while in the holster. Guide wires are double ended, stainless steel, with a 3mm J end and a floppy straight end. Product includes a vessel dilator, sheath introducer and guide wire (also available without guide wire).

## SHEATH INTRODUCERS Packaged in holster; 5 per box

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Hub Color
<b>PSI-4F-11-035</b>	4F	11cm	0.035" (0.89mm) x 50cm	red
<b>PSI-5F-11-035</b>	5F	11cm	0.035" (0.89mm) x 50cm	gray
<b>PSI-6F-11-035</b>	6F	11cm	0.035" (0.89mm) x 50cm	green
<b>PSI-7F-11-035</b>	7F	11cm	0.035" (0.89mm) x 50cm	orange
<b>PSI-8F-11-035</b>	8F	11cm	0.035" (0.89mm) x 50cm	blue
<b>PSI-5F-11-038</b>	5F	11cm	0.038" (0.97mm) x 50cm	gray
<b>PSI-6F-11-038</b>	6F	11cm	0.038" (0.97mm) x 50cm	green
<b>PSI-7F-11-038</b>	7F	11cm	0.038" (0.97mm) x 50cm	orange
<b>PSI-8F-11-038</b>	8F	11cm	0.038" (0.97mm) x 50cm	blue
<b>PSI-4F-11</b>	4F	11cm	N/A	red
<b>PSI-5F-11</b>	5F	11cm	N/A	gray
<b>PSI-6F-11</b>	6F	11cm	N/A	green
<b>PSI-7F-11</b>	7F	11cm	N/A	orange
<b>PSI-8F-11</b>	8F	11cm	N/A	blue



# Prelude® Sheath Introducers 23cm

The Prelude® Sheath Introducer is packaged in a holster that keeps all components secure, especially the guide wire. Components can be flushed while in the holster. Guide wires are double ended, stainless steel, with a floppy 3mm J end and a floppy straight end. Product includes a vessel dilator, sheath introducer and guide wire (also available without guide wire).

## SHEATH INTRODUCERS Packaged in holster; 5 per box

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Hub Color
<b>PSI-4F-23-035</b>	4F	23cm	0.035" (0.89mm) x 80cm	red
<b>PSI-5F-23-035</b>	5F	23cm	0.035" (0.89mm) x 80cm	gray
<b>PSI-6F-23-035</b>	6F	23cm	0.035" (0.89mm) x 80cm	green
<b>PSI-7F-23-035</b>	7F	23cm	0.035" (0.89mm) x 80cm	orange
<b>PSI-8F-23-035</b>	8F	23cm	0.035" (0.89mm) x 80cm	blue
<b>PSI-5F-23-038</b>	5F	23cm	0.038" (0.97mm) x 80cm	gray
<b>PSI-6F-23-038</b>	6F	23cm	0.038" (0.97mm) x 80cm	green
<b>PSI-7F-23-038</b>	7F	23cm	0.038" (0.97mm) x 80cm	orange
<b>PSI-8F-23-038</b>	8F	23cm	0.038" (0.97mm) x 80cm	blue
<b>PSI-4F-23</b>	4F	23cm	N/A	red
<b>PSI-5F-23</b>	5F	23cm	N/A	gray
<b>PSI-6F-23</b>	6F	23cm	N/A	green
<b>PSI-7F-23</b>	7F	23cm	N/A	orange
<b>PSI-8F-23</b>	8F	23cm	N/A	blue



PERCUTANEOUS ACCESS

# Prelude PRO™ Sheath Introducers & Prelude ACT™ Sheath Introducers

Prelude Pro™ Sheath tubing offers increased resilience and additional support during procedure. Prelude ACT™ (Activated Clotting Time) Sheaths allow easy blood draw through sheath side arm. Packaged in holster; 5 per box.

## PRO SHEATH INTRODUCERS

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Hub Color
<b>PRO-4F-11-035</b>	4F	11cm	0.035" (0.89mm) x 50cm	red
<b>PRO-5F-11-035</b>	5F	11cm	0.035" (0.89mm) x 50cm	gray
<b>PRO-6F-11-035</b>	6F	11cm	0.035" (0.89mm) x 50cm	green
<b>PRO-7F-11-035</b>	7F	11cm	0.035" (0.89mm) x 50cm	orange
<b>PRO-8F-11-035</b>	8F	11cm	0.035" (0.89mm) x 50cm	blue
<b>PRO-5F-11-038</b>	5F	11cm	0.038" (0.97mm) x 50cm	gray
<b>PRO-6F-11-038</b>	6F	11cm	0.038" (0.97mm) x 50cm	green
<b>PRO-7F-11-038</b>	7F	11cm	0.038" (0.97mm) x 50cm	orange
<b>PRO-8F-11-038</b>	8F	11cm	0.038" (0.97mm) x 50cm	blue
<b>PRO-4F-11</b>	4F	11cm	N/A	red
<b>PRO-5F-11</b>	5F	11cm	N/A	gray
<b>PRO-6F-11</b>	6F	11cm	N/A	green
<b>PRO-7F-11</b>	7F	11cm	N/A	orange
<b>PRO-8F-11</b>	8F	11cm	N/A	blue



## ACT SHEATH INTRODUCERS All ACT Sheaths use PRELUDE PRO™ Tubing

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Hub Color
<b>PSI-4F-11-035ACT</b>	4.5F	11cm	0.035" (0.89mm) x 50cm	red
<b>PSI-5F-11-038ACT</b>	5.5F	11cm	0.038" (0.97mm) x 50cm	gray
<b>PSI-6F-11-038ACT</b>	6.5F	11cm	0.038" (0.97mm) x 50cm	green
<b>PSI-7F-11-038ACT</b>	7.5F	11cm	0.038" (0.97mm) x 50cm	orange
<b>PSI-8F-11-038ACT</b>	8.5F	11cm	0.038" (0.97mm) x 50cm	blue
<b>PSI-4F-11ACT</b>	4.5F	11cm	N/A	red
<b>PSI-5F-11ACT</b>	5.5F	11cm	N/A	gray
<b>PSI-6F-11ACT</b>	6.5F	11cm	N/A	green
<b>PSI-7F-11ACT</b>	7.5F	11cm	N/A	orange
<b>PSI-8F-11ACT</b>	8.5F	11cm	N/A	blue

PERCUTANEOUS ACCESS

## Prelude® Sheath Introducers with Merit Advance® Needle

The Prelude® Sheath Introducer with Merit Advance needle includes the sheath introducer, vessel dilator, guide wire and needle all packaged in a holster. Packaged in holster; 5 per box.

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Needle Gauge/length	Hub Color
<b>PSI-5F-11-035-18G</b>	5F	11cm	0.035" (0.89mm) x 50cm	18G x 7cm	gray
<b>PSI-6F-11-035-18G</b>	6F	11cm	0.035" (0.89mm) x 50cm	18G x 7cm	green
<b>PSI-7F-11-035-18G</b>	7F	11cm	0.035" (0.89mm) x 50cm	18G x 7cm	orange
<b>PSI-8F-11-035-18G</b>	8F	11cm	0.035" (0.89mm) x 50cm	18G x 7cm	blue
<b>PSI-5F-11-038-18G</b>	5F	11cm	0.038" (0.97mm) x 50cm	18G x 7cm	gray
<b>PSI-6F-11-038-18G</b>	6F	11cm	0.038" (0.97mm) x 50cm	18G x 7cm	green
<b>PRO-5F-11-038-18G</b>	5F*	11cm	0.038" (0.97mm) x 50cm	18G x 7cm	gray
<b>PRO-6F-11-038-18G</b>	6F*	11cm	0.038" (0.97mm) x 50cm	18G x 7cm	green

\* Prelude PRO tubing



# Prelude® Sheath Introducers with Marker Tip

The Prelude® Sheath Introducer with Marker Tip is available in 11cm (Prelude PRO™ tubing) and 23cm (standard Prelude tubing). The sheaths have a highly radiopaque marker tip allowing for increased visibility of the sheath tip using fluoroscopy. Products are available with and without double ended, stainless steel guide wire. Packaged in holster; 5 per box.

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Hub Color
<b>PRO-4F-11-035MT</b>	4F	11cm	0.035" (0.89mm) x 50cm	red
<b>PRO-5F-11-035MT</b>	5F	11cm	0.035" (0.89mm) x 50cm	gray
<b>PRO-6F-11-035MT</b>	6F	11cm	0.035" (0.89mm) x 50cm	green
<b>PRO-7F-11-035MT</b>	7F	11cm	0.035" (0.89mm) x 50cm	orange
<b>PRO-8F-11-035MT</b>	8F	11cm	0.035" (0.89mm) x 50cm	blue
<b>PRO-5F-11-038MT</b>	5F	11cm	0.038" (0.97mm) x 50cm	gray
<b>PRO-6F-11-038MT</b>	6F	11cm	0.038" (0.97mm) x 50cm	green
<b>PRO-7F-11-038MT</b>	7F	11cm	0.038" (0.97mm) x 50cm	orange
<b>PRO-8F-11-038MT</b>	8F	11cm	0.038" (0.97mm) x 50cm	blue
<b>PRO-4F-11MT</b>	4F	11cm	N/A	red
<b>PRO-5F-11MT</b>	5F	11cm	N/A	gray
<b>PRO-6F-11MT</b>	6F	11cm	N/A	green
<b>PRO-7F-11MT</b>	7F	11cm	N/A	orange
<b>PRO-8F-11MT</b>	8F	11cm	N/A </td <td>blue</td>	blue
<b>PSI-4F-23-035MT</b>	4F	23cm	0.035" (0.89mm) x 80cm	red
<b>PSI-5F-23-035MT</b>	5F	23cm	0.035" (0.89mm) x 80cm	gray
<b>PSI-6F-23-035MT</b>	6F	23cm	0.035" (0.89mm) x 80cm	green
<b>PSI-7F-23-035MT</b>	7F	23cm	0.035" (0.89mm) x 80cm	orange
<b>PSI-8F-23-035MT</b>	8F	23cm	0.035" (0.89mm) x 80cm	blue
<b>PSI-5F-23-038MT</b>	5F	23cm	0.038" (0.97mm) x 80cm	gray
<b>PSI-6F-23-038MT</b>	6F	23cm	0.038" (0.97mm) x 80cm	green
<b>PSI-7F-23-038MT</b>	7F	23cm	0.038" (0.97mm) x 80cm	orange
<b>PSI-8F-23-038MT</b>	8F	23cm	0.038" (0.97mm) x 80cm	blue
<b>PSI-4F-23MT</b>	4F	23cm	N/A	red
<b>PSI-5F-23MT</b>	5F	23cm	N/A	gray
<b>PSI-6F-23MT</b>	6F	23cm	N/A	green
<b>PSI-7F-23MT</b>	7F	23cm	N/A	orange
<b>PSI-8F-23MT</b>	8F	23cm	N/A	blue



PERCUTANEOUS ACCESS

## Prelude® Sheath Introducers - Mini Access Kits

The Prelude® Sheath Introducers - Mini Access Kits provide access to the vascular system to facilitate the placement of guide wires. The sheaths are available in three lengths: 7cm, 11cm and 23cm; and are packaged in a holster with a sheath, vessel dilator, guide wire and needle. Products can be used for pediatrics, transradial, brachial and femoral access. Packaged in holster; 5 per box.

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Tip	Hub Color
<b>PSI-4F-7-018</b>	4F	7cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	red
<b>PSI-5F-7-018</b>	5F	7cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	gray
<b>PSI-6F-7-018</b>	6F	7cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	green
<b>PSI-4F-7-018NT4</b>	4F	7cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	red
<b>PSI-5F-7-018NT4</b>	5F	7cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	gray
<b>PSI-6F-7-018NT4</b>	6F	7cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	green
<b>PSI-4F-11-018</b>	4F	11cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	red
<b>PSI-5F-11-018</b>	5F	11cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	gray
<b>PSI-6F-11-018</b>	6F	11cm	0.018" (0.46mm) x 40cm	SS/SS	21G x 4cm	green
<b>PSI-4F-11-018NT4</b>	4F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	red
<b>PSI-5F-11-018NT4</b>	5F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	gray
<b>PSI-6F-11-018NT4</b>	6F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 4cm	green
<b>PRO-4F-11-018NT7*</b>	4F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	red
<b>PRO-5F-11-018NT7*</b>	5F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	gray
<b>PRO-6F-11-018NT7*</b>	6F	11cm	0.018" (0.46mm) x 40cm	NT/PT	21G x 7cm	green
<b>PSI-4F-23-018NT4</b>	4F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 4cm	red
<b>PSI-5F-23-018NT4</b>	5F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 4cm	gray
<b>PSI-6F-23-018NT4</b>	6F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 4cm	green
<b>PSI-4F-23-018NT7</b>	4F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 7cm	red
<b>PSI-5F-23-018NT7</b>	5F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 7cm	gray
<b>PSI-6F-23-018NT7</b>	6F	23cm	0.018" (0.46mm) x 80cm	NT/PT	21G x 7cm	green

SS = Stainless Steel PT = Platinum Tip NT = Nitinol

Guide wire is stainless steel; double ended (3mmJ and floppy).

\* Catalog Numbers that begin with PRO use Prelude PRO tubing

<b>PSI-4F-7-025</b>	4F	7cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	red
<b>PSI-5F-7-025</b>	5F	7cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	gray
<b>PSI-6F-7-025</b>	6F	7cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	green
<b>PSI-4F-11-025</b>	4F	11cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	red
<b>PSI-5F-11-025</b>	5F	11cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	gray
<b>PSI-6F-11-025</b>	6F	11cm	0.025" (0.64mm) x 50cm	-	20G x 4cm	green



# Prelude® Short Sheath Introducers

Prelude® Short Sheath Introducer provides access to dialysis fistulas and grafts. Can be used for declot, angioplasty and thrombolysis procedures, as well as temporary hemodialysis. Will accommodate flow rates up to 600ml/min. Includes removable stopcock. Packaged 5 per box.

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Needle Gauge/Length	Hub Color
<b>PSS-5F-4-018MT</b>	5F	4cm	0.018" (0.46mm) x 40cm*	21G x 4cm	gray
<b>PSS-6F-4-018MT</b>	6F	4cm	0.018" (0.46mm) x 40cm	21G x 4cm	green
<b>PSS-7F-4-018MT</b>	7F	4cm	0.018" (0.46mm) x 40cm	21G x 4cm	orange
<b>PSS-8F-4-018MT</b>	8F	4cm	0.018" (0.46mm) x 40cm	21G x 4cm	blue
<b>PSS-5F-4-038MT</b>	5F	4cm	0.038" (0.97mm) x 50cm**	–	gray
<b>PSS-6F-4-038MT</b>	6F	4cm	0.038" (0.97mm) x 50cm	–	green
<b>PSS-7F-4-038MT</b>	7F	4cm	0.038" (0.97mm) x 50cm	–	orange
<b>PSS-8F-4-038MT</b>	8F	4cm	0.038" (0.97mm) x 50cm	–	blue
<b>PSS-5F-4MT</b>	5F	4cm	N/A	–	gray
<b>PSS-6F-4MT</b>	6F	4cm	N/A	–	green
<b>PSS-7F-4MT</b>	7F	4cm	N/A	–	orange
<b>PSS-8F-4MT</b>	8F	4cm	N/A	–	blue

\* 0.018" Guide wire is stainless steel with straight stainless steel flexible tip.

\*\* 0.038" Guide wire is stainless steel; double ended (3mm J and floppy).

