



Marquis[®] Stopcocks

- 254 | 1-Way Stopcocks
- 255 | 3-Way Stopcocks
- 257 | 4-Way Stopcocks
- 259 | Stopcocks with Extension Tubing

SEE RELATED
PRODUCTS



Tubing
223



Custom Kits
87



Caps & Covers
103

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

1-Way Stopcocks

Handle: Acetal Material: Polycarbonate

1200 PSI



Fixed Male Luer

U1FLP



Swivel Nut

U1SNP



Right Rotating Male Luer

U1RP



Slip Luer

U1SLP



Fixed Male Luer

H1FLP



Swivel Nut

H1SNP



Right Rotating Male Luer

H1RP

1050 PSI



Fixed Male Luer

M1FLP



Swivel Nut

M1SNP



Slip Luer

M1SLP

500 PSI



Fixed Male Luer

S1FLP



Swivel Nut

S1SNP



Slip Luer

S1SLP

200 PSI

3-Way Stopcocks

Handle: Acetal Material: Polycarbonate

1200 PSI



Slip Luer
U3SLP



Fixed Male Luer
U3FLP



Right Rotating Male Luer
U3RRP



Swivel Nut
U3SNP

1050 PSI



Slip Luer
H3SLP



Fixed Male Luer
H3FLP



Left Rotating Male Luer
H3RLP



Right Rotating Male Luer
H3RRP



Swivel Nut
H3SNP

500 PSI



Slip Luer
M3SLP



Fixed Male Luer
M3FLP



Left Rotating Male Luer
M3RLP



Right Rotating Male Luer
M3RRP



Swivel Nut
M3SNP

3-Way Stopcocks

Handle: Acetal Material: Polycarbonate



Swivel Nut
S3SNP



Slip Luer
S3SLP



All Female Luer
S3FP



Fixed Male Luer
S3FLP

200 PSI



Swivel Nut, Large Bore
L3SNP

50 PSI

4-Way Stopcocks

Handle: Acetal Material: Polycarbonate

1200 PSI



Right Rotating
Male Luer
U4RRP

Swivel Nut
U4SNP

Fixed Male Luer
U4FLP

Slip Luer
U4SLP

1050 PSI



Fixed Male Luer
H4FLP

Right Rotating Male Luer
H4RRP

Swivel Nut
H4SNP

500 PSI



Right Rotating
Male Luer
M4RRP

Swivel Nut
M4SNP

Fixed Male Luer
M4FLP

Right Rotating Male
Luer w/ On Handle
MN4RRP

4-Way Stopcocks

Handle: Acetal Material: Polycarbonate



50 PSI, Swivel Nut,
Large Bore

L4SNP



200 PSI,
Swivel Nut

S4SNP



200 PSI,
Fixed Male Luer

S4FLP



200 PSI,
Slip Luer

S4SLP

50 to 200 PSI