

## Stopcocks with Extension Tubing

Merit offers a variety of stopcocks with tubing extensions to meet various project needs. In addition to the wide variety listed below, we can also offer a solution tailored to your needs. Contact your sales representative or OEM customer service for more information.

### Stopcocks with Extension Tubing Examples:



3-way Stopcock  
Standard Pressure  
0.110" Tubing OD  
DEHP-free PVC Tubing  
5.625" Extension

**501857001**



3-way Stopcock  
Standard Pressure  
0.158" Tubing OD  
DEHP-free PVC Tubing  
6.88" Extension

**501869001**



3-way Stopcock  
Standard Pressure  
0.140" Tubing OD  
TPU Tubing  
8" Extension

**502226001**



1-way Stopcock  
Medium Pressure  
0.125 Tubing OD  
PVC Tubing  
6.25" Extension with  
Male Luer

**K04-01006P**



1-way Stopcock  
Medium Pressure  
0.125 Tubing OD  
PVC Tubing  
6.25" Extension  
with Male Luer

**K04-01092P**



1-way Stopcock  
Standard Pressure  
0.215" Tubing OD  
PVC Tubing  
4" Extension

**K12-05624**