

Tear Away/ClassicSheath™ Introducers



Sheath ID: 0.083"
 Size: 6F, Color: Green
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1006



Sheath ID: 0.095"
 Size: 7F, Color: Orange
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1007



Sheath ID: 0.105"
 Size: 8F, Color: Blue
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1008



Sheath ID: 0.112"
 Size: 8.5F, Color: Light Blue
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-10008.5



Sheath ID: 0.118"
 Size: 9F, Color: Black
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1009



Sheath ID: 0.118"
 Size: 9.5F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1009.5

Tear Away/ClassicSheath™ Introducers



Sheath ID: 0.131"
 Size: 10F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1010



Sheath ID: 0.138"
 Size: 10.5F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1010.5



Sheath ID: 0.144"
 Size: 11F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1011



Sheath ID: 0.158"
 Size: 12F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1012



Sheath ID: 0.164"
 Size: 12.5F, Color: White
 Length: 13cm
 Guide Wire Diameter/Length:
 .038" x 45cm

CLS-1012.5



Sheath ID: 0.092"
 Size: 7F, Color: Orange
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2507

PERCUTANEOUS
 ACCESS

Tear Away/ClassicSheath™ Introducers



Sheath ID: 0.105"
 Size: 8F, Color: Blue
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2508



Sheath ID: 0.118"
 Size: 9F, Color: Black
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2509



Sheath ID: 0.131"
 Size: 10F, Color: White
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2510



Sheath ID: 0.138"
 Size: 10.5F, Color: White
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2510.5



Sheath ID: 0.144"
 Size: 11F, Color: White
 Length: 25cm
 Guide Wire Diameter/Length:
 .038" x 80cm

CLS-2511