

Transducers

- Reliable Precision
- Clear Fluid Path
- Customized Configurations also Available

SPECIFICATIONS

| | |
|--------------------------|----------------------------|
| Operating Pressure Range | -50 to +300 mmHg |
| Sensitivity | 5 μ V/V/mmHg, \pm 1% |
| Zero Offset | 0 \pm 25 mmHg |

Single use. Sold packaged non-sterile. For sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED
PRODUCTS



Syringes
195



Stopcocks
253



Manifolds
249



MERIT SENSOR

Merit Medical Systems' wholly owned subsidiary, Merit Sensor, designs and manufactures pressure sensors included in the Merit Transducer as well as Merit's digital inflation devices. See page 157 or go to www.MeritSensor.com

Meritrans® Blood Pressure Transducers



With Non-Winged Housing and 3-Way and 1-Way Stopcocks with 48" (122cm) cable

K11-00202P



With Non-Winged Housing and 3-Way Stopcock (Med Pressure, Rotating Male Adapter) and 1-Way Stopcock with 48" (122cm) cable

K11-00392AP



With Winged Housing and 2 Bonded Standard Pressure 1-Way Stopcocks with 12" (30cm) cable

K11-00098AP



With Non-Winged Housing and 3-Way Stopcock bonded to both ends with 48" (122cm) cable

K11-00144AP



With Non-Winged Housing and Bonded Standard Pressure 1-Way Stopcock with 48" (122cm) cable

K11-00137AP



With Winged Housing and Bonded Standard Pressure 3-Way and 1-Way Stopcocks with 12" (30cm) cable

K11-00100AP

SPECIALTY
DEVICES

Meritrans® Blood Pressure Transducers



With Non-Winged Housing with
12" (30cm) cable

MER212P



With Non-Winged Housing with
48" (122cm) cable

MER200P



With Winged Housing with
12" (30cm) cable

MER100P



With Winged Housing and 2 Bonded
Standard Pressure 3-Way Stopcocks with
12" (30cm) cable

K11-00099AP



With Winged Housing and Bonded
Standard Pressure 1-Way and 3-Way
Stopcocks with 12" (30cm) cable

K11-00096AP



With Winged Housing with
48" (122cm) cable

MER148P

SPECIALTY
DEVICES

Meritrans® Blood Pressure Transducers



Medium Pressure Transducer
0 - 100 PSI
501669001



High Pressure Transducer
0 - 440 PSI
501750001

Transducers Flush Device



Adult 3cc/hr flow rate.
White
270428001



Pediatric 30cc/hr flow rate.
Yellow
270428002