

Merit X-Change™ Syringes

Plunger Material: ABS
 Material: Polycarbonate
 Tip Material: Silicone



20ml X-Change™, clear, with male luer tip
EX20-SMP



20ml X-Change™, gamma, with male luer tip
EX20G-SMP



20ml X-Change™, clear, with dead-end cap
EX20-SDP



20ml X-Change™, gamma, with dead-end cap
EX20G-SDP

Also Available Separately

Male luer tips



Clear
702395001



Gamma
702395003

Dead-end caps



Clear
702395002



Gamma
702395004

CUSTOMIZING

From custom pad printing on the barrel to custom colored plungers or custom barrel tip configurations, Merit Medical can work with your product design team to customize the X-Change™ to fit your needs.

