




Confidential Information. Property of Merit Medical Systems, Inc.
Reproduction or use without written consent is prohibited.

DRAWING IS FOR REFERENCE ONLY



MALE LUER LOCK FITTING PER ISO 594



NOTES:

MATERIALS:

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (WHITE)
- SYRINGE TIP: SILICONE (BLACK)

	MERITMEDICAL™	OEM DIVISION
1600 West Merit Parkway • South Jordan, Utah 84095		
PART NUMBER: CCS601		
TITLE: Syringe, CCS, 12mL, Ring Grip, Ring Plunger		
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY.		
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.		