




**NOTES:**  
**MATERIALS:**

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (WHITE)
- SYRINGE TIP: SILICONE (BLACK)
- ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS  
(CONTACT A BUSINESS MANAGER)

<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

DRAWING IS FOR REFERENCE ONLY



	<b>OEM</b>	1600 West Merit Parkway South Jordan, Utah 84095
DRAWING TITLE: SYRINGE, CCS, 12ML, FIXED MALE LUER, PALM PAD PLUNGER		
PART NUMBER:	CCSB701P	REV.: N/A
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.		DATE: 11/03/2014