


**NOTES:**

**MATERIALS:**

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (YELLOW)
- SYRINGE TIP: SILICONE (BLACK)
- ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS  
(CONTACT A BUSINESS MANAGER)  
<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

	OEM 1600 West Merit Parkway South Jordan, Utah 84095
DRAWING TITLE: SYRINGE, MEDALLION, 20 mL, WING GRIP, YELLOW	PART NUMBER: MSSW21-YP REV.: N/A
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.	DATE: 12/12/2016