



STD LUER TAPER
PER ISO 594
LUER TAPER PER
ISO 80369-7 SEE NOTES



NOTES:

MATERIAL: POLYPROPYLENE (YELLOW)
ISO 80369-7 IN PROCESS: CONTACT
CUSTOMER SERVICE

<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

MERITMEDICAL OEM DIVISION
1600 West Merit Parkway • South Jordan, Utah 84095

DRAWING TITLE: COVER, FEMALE, NON-VENTED, YELLOW

PART NUMBER: 101020004

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY.
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.