


**NOTES:**

**MATERIALS:**

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (YELLOW)
- SYRINGE TIP: SILICONE (BLACK)
- ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS (CONTACT A BUSINESS MANAGER)  
<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

 <b>MERITMEDICAL</b>	<b>OEM</b>	1600 West Merit Parkway South Jordan, Utah 84095
DRAWING TITLE: SYRINGE, MEDALLION, 20ML, YELLOW		
PART NUMBER: MSS121-YP	REV.: N/A	DATE: 01/23/2015
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.		