



# NOTES:

## MATERIALS:

- BARREL: POLYCARBONATE (CLEAR)
  - PLUNGER: ABS (YELLOW)
  - SYRINGE TIP: SILICONE (BLACK)
  - ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS (CONTACT A BUSINESS MANAGER)
- <https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

|   |                      |                     |
|---|----------------------|---------------------|
|                                  | <b>MERITMEDICAL™</b> | <b>OEM DIVISION</b> |
| 1600 West Merit Parkway • South Jordan, Utah 84095  |                      |                     |
| DRAWING TITLE: SYRINGE, MEDALLION, 10 mL, YELLOW  |                      |                     |
| PART NUMBER: MSS111-YP  |                      |                     |
| UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY.<br>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. |                      |                     |