


ISO 594  
THREAD  
LUER THREAD PER ISO  
80369-7; SEE NOTES



SECTION A-A

NOTES

MATERIAL: ABS, YELLOW  
ISO 80369-7 IN PROCESS: CONTACT  
CUSTOMER SERVICE

<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

	<b>MERITMEDICAL™</b>	<b>OEM DIVISION</b>
<small>1600 West Merit Parkway • South Jordan, Utah 84095</small>		
DRAWING TITLE: CAP, MALE, NON-VENTED, YELLOW		
PART NUMBER: 101041004		
<small>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.</small>		