


**NOTES:**

**MATERIALS:**


- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (PURPLE)
- SYRINGE TIP: SILICONE (BLACK)
- ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS (CONTACT A BUSINESS MANAGER)  
<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>


DRAWING IS FOR REFERENCE ONLY



MALE LUER LOCK FITTING  
PER ISO 594-2

MALE LUER LOCK FITTING  
PER ISO 80369-7; 2021 SEE NOTE



	<b>MERITMEDICAL</b>	<b>OEM</b>	1600 West Merit Parkway South Jordan, Utah 84095
DRAWING TITLE: SYRINGE, MEDALLION, 1ML, PURPLE			
PART NUMBER: MSS011-PRP		REV.: N/A	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.			DATE: 01/26/2015