


NOTES:

MATERIALS:

- BARREL: POLYCARBONATE (CLEAR)
 - PLUNGER: ABS (YELLOW)
 - SYRINGE TIP: SILICONE (BLACK)
 - ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS (CONTACT A BUSINESS MANAGER)
- <https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

		OEM		1600 West Merit Parkway South Jordan, Utah 84095	
DRAWING TITLE: SYRINGE, MEDALLION, 10 mL, WING GRIP, YELLOW					
PART NUMBER:		MSSW11-YP			REV.: N/A
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.				DATE:	1/30/2017