



NOTES:
MATERIALS:

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (PURPLE)
- SYRINGE TIP: SILICONE (BLACK)
- ISO 80369-7 COMPLIANCE TRANSITION IN PROCESS (CONTACT YOUR BUSINESS MANAGER)

<https://meritoem.com/about/iso-80369-small-bore-connector-standard-iso-80369-7-for-vascular-and-intravenous-applications/>

	MERITMEDICAL™	OEM DIVISION
<small>1600 West Merit Parkway • South Jordan, Utah 84095</small>		
DRAWING TITLE: SYRINGE, MEDALLION, 10 mL, PURPLE, LIDOCAINE		
PART NUMBER: K01-04086P		
<small>UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY. UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.</small>		