



**NOTES:**

**MATERIALS:**

- BARREL: POLYCARBONATE (CLEAR)
- PLUNGER: ABS (YELLOW)
- SYRINGE TIP: SILICONE (BLACK)



**MERITMEDICAL™**

DRAWING TITLE: 6 mL, YELLOW, MEDALLION SLIP LUER SYRINGE

PART NUMBER: K01-02044P

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE FOR REFERENCE ONLY.  
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.