

VacLok AT™

Barrel Material: Polycarbonate Tip Material: Silicone
Plunger Material: Acrylic Copolymer

A History of Innovation, A Future of Advancement

Designed to hold a vacuum without user interface, Merit's VacLok AT™ locks in any of the 90 positions along the length of the barrel – holding negative pressure with ease and reducing hand fatigue, making it the perfect syringe for creating negative vacuum pressure.

SEE THE DIFFERENCE

- Accurate filling
- Clear Polycarbonate barrel
- Large, stable finger wings
- Bold graduations
- Silicone plunger tip
- Dual use – can be used as a vacuum syringe or as a standard syringe

FEEL THE DIFFERENCE

- Smooth glass plunger action
- Even, steady operation
- Larger grip
- Skid resistant thumb pad
- Plunger design holds negative pressure with ease and won't inadvertently disengage



Step 1: To activate vacuum pressure, pull plunger back to desired location. Twist plunger clockwise until it stops and locks into place.



Step 2: To disconnect vacuum pressure, twist plunger counter-clockwise until plunger stops.



20ml, 0-90 positions

VAC220