

Cultura™

FLEX SWAB

Optimal Collection. Improved Comfort.

Introducing the new Cultura™
Flex Swab from Merit Medical.
Developed as a critical part of
streamlining the process of
collecting specimen samples
for a nasopharyngeal swab.

Our propriety Flex Swab tip is designed with a Nylon fiber flocking
for patient comfort and provides superior collection and release
of specimens. The Cultura™ Flex Swab can be custom labeled
and package as either bulk non-sterile or sterile. Manufactured at
Merit's Global Headquarters in South Jordan, Utah.



Cultura™ FLEX SWAB



FEATURES & BENEFITS


- Our propriety Flex Swab tip is designed with a Nylon fiber flocking for patient comfort and provides superior collection and release of specimens.
- The flex swab is designed for use with our twist on cap, the swab can be secured and locks into place during transport, reducing risk of contamination when removing the swab.
- Individually packaged as a sterile unit. Other options include custom and private labeling.
- Features two breakpoint option lengths of 80 mm and 100 mm from the swab tip.
- Intended use as a nasopharyngeal swab.

SPECIFICATIONS


- Handle/Shaft molded with high-impact polystyrene.
- Swab tip is 16 mm in length and 3.3 mm wide.
- Flex Swap tip is Nylon fiber flocked.
- FDA 510K and CE registered.

Catalog Number	Description
C100F	100 mm breakpoint, 150 mm overall length
C80F	80 mm breakpoint, 150 mm overall length

Carton Units	Shipper Units
75 units per carton	300 units per shipper
250 units per carton	1000 units per shipper



Capture-cap aligns swab and allows the swab to self-seat with minimal insertion force



Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.

MERIT MEDICAL | OEM
COMPONENTS | DEVICES | POSSIBILITIES™

Merit Medical OEM
1600 West Merit Parkway
South Jordan, Utah 84095 USA

To place an order in the
US, contact:
Toll Free +1 800 637 4839
Tel. +1 801 208 4313
Fax +1 801 253 6988
Email Oemorder@Merit.com

To place an order in
Europe and Asia, contact:
Tel. +31 43 358 82 24
(Maastricht, The Netherlands)
Fax +31 43 201 09 06
Email Oemorder@Merit.com

MeritOEM.com