

Catheter Hubs

Material: Polycarbonate



4F Hub, White
100941200



5F Hub, White
100941205



6F Hub, White
.050" thru-hole
100941211



6F Hub, White
100941212



8F Hub, White
100941218



4F Hub, Clear
100942220



5F Hub, Clear
100942221



6F Hub, Clear,
.083" Inner Lumen
100942222




7F Hub, Clear
100942223



7F Hub, Clear
100942224

Hub details. Measurements in inches.



Catalog Number	A	B	C	D
100941200	0.087	0.090	0.058	0.056
100941205	0.111	0.090	0.071	0.069
100941211	0.120	0.125	0.083	0.050
100941212	0.012	0.0125	0.083	0.052
100941218	0.135	0.125	0.109	0.107

C $^{+.002}_{-.001}$ \varnothing at entrance straightline tapers to D $^{+.001}_{-.001}$ \varnothing at bottom



Custom colors
and printing
options available.



5F Hub, White
0.045 inner lumen
Printed

590012039



Reverse



6F Hub, White
0.050 inner lumen
Printed

590012057



Reverse